

PDI®



Material Safety Data Sheet

PDI® - PROFESSIONAL DISPOSABLES INTERNATIONAL
 THE HEALTHCARE DIVISION OF NICE-PAK PRODUCTS, INC
 TWO NICE-PAK PARK
 ORANGETBURG NY 10962-1376
 845-365-1700
 REVISION DATE: 3/03

I. PRODUCT IDENTIFICATION

PRODUCT/TRADE NAME: **SANI-DEX® ALC Antimicrobial Alcohol Hand Wipes – Item# D20300**

HAZARD RATING(NFPA)

HEALTH: 1
 FLAMMABILITY: 3
 REACTIVITY: 0
 SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F daytime) at other times, contact the local Poison Control Center.

CHEMICAL NAME: Proprietary Mixture

II HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200

Hazardous Ingredients	% w/w	ACGIH TLV	CAS NUMBER
SD ALCOHOL 40 (190 Proof)	65.3	1000 PPM	64-17-5

III-PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: Towelette saturated with a clear, colorless liquid with an alcohol odor

BOILING POINT: N/A

FLASH POINT: 74⁰ F

VAPOR DENSITY: N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: Complete

SPECIFIC GRAVITY (H₂O = 1): 0.890

IV FIRE & EXPLOSION HAZARD DATA

FLASH POINT (Method Used): 74⁰ F Tag Closed Cup : LEL: N/A UEL: N/A

EXTINGUISHING MEDIA: Carbon dioxide or alcohol type foam

SPECIAL FIRE FIGHTING PROCEDURES: Handle as flammable liquid. Use respiratory protection for fire fighting personnel.

UNUSUAL FIRE AND EXPLOSION HAZARDS: NFPA Class 3 Flammability

V-REACTIVITY DATA

STABILITY: Stable

CONDITIONS TO AVOID: None Known

INCOMPATIBILITY: None

HAZARDOUS DECOMPOSITION OR BYPRODUCT: None Known

POLYMERIZATION: Will not occur

CONDITIONS TO AVOID: None known

PDI[®]



SANI-DEX[®] ALC Item #D20300

VI-HEALTH HAZARD DATA

EFFECTS OF OVER EXPOSURE:

SKIN Formulated for topical application
EYES : Splashes may cause irritation
INHALATION: None known
INGESTION Not a normal route of entry

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: Discontinue use if rash
EYE CONTACT: Flush with cool water for 15 minutes
INHALATION: None
INGESTION: Consult Physician
TARGET ORGANS: None known

VII-SPILL AND DISPOSAL PROCEDURE

SPILL CONTROL: Pick up with absorbent media
WASTE DISPOSAL METHOD: Obey local, state and federal regulations
HANDLING AND STORAGE: Store away from ignition sources such as heat, sparks, open flame, etc...

VIII-CONTROL MEASURES/PROTECTION

RESPIRATION: None required
VENTILATION: Recommended
PROTECTIVE GLOVES: None required
EYE PROTECTION: Only if splashing is expected
HYGIENIC PRACTICES: Good housekeeping practices should be followed
OTHER: None

IX TRANSPORT/SHIPPING

DOT SHIPPING NAME: CONSUMER COMMODITY
TECHNICAL SHIPPING NAME: Not applicable
DOT SHIPPING CLASSIFICATION: DOT not regulated (except by boat)
DOT ID NO.: Not applicable
DOT LABEL REQUIREMENTS: Not applicable
UN/NA NUMBER REGULATIONS: Not applicable
REPORTABLE QUANTITY: Not applicable

X DISCLAIMER

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.

N/A = NOT AVAILABLE